

Metron PC System Requirements – Minimum

Below are the minimum requirements for the MetronDVM software to operate but may not display images at the speed and image quality you might need.

Operating System	Windows XP Pro 32-bit or Windows 7 32-bit
Graphics System and Monitor	Minimum requirement of 1024x768 pixel display with fonts at “normal” 96 DPI. Strongly recommend 1280x1024 or greater.
RAM	Minimum 1 GB for Windows XP Pro and minimum 2GB for Windows 7
Disk Space	To install Metron, the computer must have at least 500 MB of available disk space. Additional room will be needed for the images - assume needing at least 2 GB for every 100 images.
Processor Speed	Metron will typically perform better on newer/faster computers. Computers that are older have slower processor speeds. To help maximize the speed/efficiency of the Metron program and your overall image processing/viewing experience, use the fastest computer possible.
Internet Connection	High speed Internet connection availability is required on any computer that will be running the Metron product. The internet connection is critical for training & support. Also, with Metron, you can send diagnostic images electronically to a remote radiologist; this is not possible without an Internet connection.

Metron PC System Requirements - Suggested Configuration

Below are the suggested performance requirements for the MetronDVM software.

Operating System	Windows XP Pro 32-bit or Windows 7 32-bit
Graphics System and Monitor	Suggested monitor selection is at least 20 inches wide, max resolution of 1600 x 1200 or greater / 60 Hz, image contrast ratio – 800:1, and response time of 16 ms. Video card should have separate memory of 512 MB.
RAM	Computers running Windows XP Pro will perform best with at least 2 GB of RAM. Computers running Windows 7 will perform best with at least 4 GB of RAM.
Disk Space	Consider 160 - 320 GB disk as the minimum to avoid needing a larger disk later.
Processor Speed	Metron will typically perform better on newer/faster computers. Computers that are older have slower processor speeds. To help maximize the speed/efficiency of the Metron program and your overall image processing/viewing experience, use the fastest computer possible.
Internet Connection	High speed Internet connection availability is required on any computer that will be running the Metron product. The internet connection is critical for training & support. Also, with Metron, you can send diagnostic images electronically to a remote radiologist; this is not possible without an Internet connection.